

The **Training Checklist** serves as a guide during the training process. It is imperative that all skills set forth below are learned and understood so that proper use and implementation of the Nisus NPWT pump will be achieved.

Instructor is required to place his/her initials in the appropriate column designated below to indicate if the trainee skill was observed, or if the trainee skill was verbally explained back to the instructor.

Staff Member Name:

Date:

Home Screen		
Skill:	Skill observed:	Skill verbally explained:
Powering on pump		
Understanding date and time		
Understanding the battery indicator		
Understanding of the pressure setting		
Understanding of the lock symbol		
Understanding of the run time		
Understanding of the mode of operation read out		
Basic Features		
Skill:	Skill observed:	Skill verbally explained:
Identifying the power button		
Identifying the up arrow, down arrow, right arrow, and left arrow buttons		
Identifying the Menu/Select button		
Identifying the Exit button		
Identifying the charger port		
Identifying the USB port		
Identifying the serial number		
Attaching the battery charger		
Identifying the charging indicator light		
Therapy Menu		
Skill:	Skill observed:	Skill verbally explained:
Setting to Continuous Mode of operation		
Setting pressure for Continuous Mode of therapy		
Understanding the current mm/hg readout		
Setting to Intermittent Mode of operation		
Setting up pressure		
Setting up time		
Setting of down pressure		
Setting down time		
Understanding the current mm/hg readout		

*Forms are to be used as a guide *only*; please consult and follow applicable policies and procedures established by your company.

Alarm Menu		
Skill:	Skill observed:	Skill verbally explained:
How to turn Leak Alarm on and off		
How to turn Blockage Alarm on and off		
How to turn on and off the Troubleshoot Menu		
Understanding of the Display Log		
Setting Menu		
Skill:	Skill observed:	Skill verbally explained:
How to change language		
How to change date and time		
How to access Factory Admin. screen		
In the Factory Admin. screen – how to reset the Alert Log		
How to reset the run time		
Explain the software version		
Troubleshooting Menu		
Skill:	Skill observed:	Skill verbally explained:
Explanation of Leak troubleshooting		
Explanation of Blockage troubleshooting		
Explanation of Canister Full troubleshooting		
Canisters		
Skill:	Skill observed:	Skill verbally explained:
Explain the two different sizes of canisters		
How to attach a canister		
How to remove a canister		
Explain how often to change a canister		
Advanced Features		
Skill:	Skill observed:	Skill verbally explained:
How to lock and unlock pump		
Explanation of Sleep Mode		
How to turn Sleep Mode on and off		
Explanation of what alarms can be muted		
How to mute alarms		
Explanation of Canister Full alarm		
Explanation of the two different battery alarms		

I have reviewed this document in its entirety and by signing below, attest that I understand how to properly use and operate the Nisus NPWT Pump.

Signature of staff member: _____ Date: _____

Signature of instructor: _____ Date: _____