



**Nisus Li-Ion Battery, 14.8volt/ 1300mah**

**Part Number: CK-2301**

**USE**

The Nisus NPWT Pump is designed to run on a lithium-ion rechargeable battery. The pump can also operate with the provided battery charger plugged in.

**INDICATIONS**

The Nisus Li-Ion Battery is indicated for:

- Use with the Nisus NPWT Pump.

**CONTRAINDICATIONS**

The Nisus Li-Ion Battery is contraindicated for:

- Not for use with any other device, medical or personal.

**WARNINGS/CAUTIONS**


- Read and understand all warnings, cautions, and directions completely and carefully before use. Failure to do so may result in serious injury.
- Do not use if the individual package is opened or if the package or the product is damaged. Carefully inspect the product for damage or defects prior to use.
- Disposal for Li-Ion Battery Pack:  
The rechargeable battery contains the following components that must be disposed of in accordance with local regulations: Lithium Ion cells.  
Caution: When the battery is worn out, insulate the connector terminals with adhesive tape or similar material before disposal.  
Misuse or improper disposal of the battery pack may cause the battery to become very hot, ignite, or explode. When disposing of a battery pack, contact your local waste disposal service provider regarding local restrictions on the disposal and recycling of batteries.

**DIRECTIONS FOR USE**

1. Secure battery into designated location inside of Nisus unit using supplied double stick adhesive strips.
2. Insert battery connector into designated circuit board location.
3. Close unit, secure with fasteners.
4. Power on device to test connection.

**STORAGE AND TRANSPORT CONDITIONS**

Store at room temperature. Avoid excessive heat and humidity.

 **Manufactured For:**  
Cork Medical Products  
8000 Castleway Drive  
Indianapolis, IN 46250 USA