



## Elite Balance Air Mattress

Nineteen individual air cells provide a unique design which keeps the lower layer of air cells constantly inflated while alternating cells inflate and deflate the upper layer.

### FEATURES OF THE ELITE BALANCE AIR

- ▶ Digital pump with light-up display allows for effortless user experience and control
- ▶ Pressure control settings allow patient to set their personalized comfort setting depending on weight
- ▶ Auto-Detection setting automatically detects the patient weight and sets the appropriate pressure output for optimum patient comfort and ease
- ▶ Pump features audible and visual low pressure, power failure and machine malfunction alarms
- ▶ Split cell-on-cell design allows the patient to remain in an upright position and keeps mattress inflated up to 24 hours if loss of power should occur
- ▶ In the event of a medical emergency, mattress features a rapid deflation CPR valve
- ▶ Standard and HD Mattresses are available with or without bolsters

HCPCS Code: **E0277**

### COMPONENT SPECS

- 19 air cells
- Cell Material: Nylon TPU
- Cover Material: Stretched Polyester with Polyurethane, replacement covers available
- 1-Year Warranty

### AVAILABLE IN 3 SIZES

Model #	Model Name	Dimensions	Max. Wt
6406-05 or 6406-05-Bolster	EBA Standard	78" (L) x 35" (W) x 8" (H)	350 lbs.
6406-08 or 6406-08-Bolster	EBA HD	78" (L) x 42" (W) x 8" (H)	550 lbs.
6406-11	EBA MAX	78" (L) x 48" (W) x 8" (H)	750 lbs.

