




Range: 10.0 to 200.0 mmHg Average Pressure: 22.4 mmHG Peak Pressure: 35.9 mmHG Area: 359.38



35" Elite Balance Air

42" Elite Balance Air HD

48" Elite Balance Air Max

Model	6406-05	6406-08	6406-11
Mattress Size	78" (L) x 35" (W) x 8" (H)	78" (L) x 42" (W) x 8" (H)	78" (L) x 48" (W) x 8" (H)
Cell Height	8"	8"	8"
Number of Cells	19	19	19
Cell Material	Nylon-PU	Nylon-PU	Nylon-TPU
Maximum weight capacity	350 lbs.	550 lbs.	750 lbs.
Pillow Section	Yes	Yes	Yes
Bolsters	Yes, 6406-05-Bolster	Yes, 6406-08-Bolster	No
Cell-on-Cell Design	Yes	Yes	Yes
Mattress Base Material	Polyester-PU	Polyester-PU	Polyester-PU
Number of Knitting Cell Retainers	19	19	19
Independent Low Air Loss	Yes	Yes	Yes
CPR	Yes	Yes	Yes
Alternating Type	A-B Alternating		
Pump Output	9.0L		18.0L
Pressure Control Type	Digital (Digital Keypad + LED display)		
Alarms	LED with Audible Alarms (Power, Low Pressure, Service)		
Pressure Ranges	25-60mmHg, ±5mmHg		
Auto Static-to-Alternate Fallback	Yes	Yes	Yes
Auto-Firm	Yes	Yes	Yes
Auto-Weight Detection	Yes	Yes	Yes
Alternating Cycle Time	Yes (10-15-20-25-30 Minute Intervals)		
Water Resistant (IP21)	Yes	Yes	Yes
HCCPS Codes	E0277	E0277	E0277
Warranty	1 Year	1 Year	1 Year