

Prevention Cushion





Treatment Cushion



Freedom Ease

Forever Plus

	Freedom Ease			Forever Plus		
	Model	Dimensions	Cells	Model	Dimensions	Cells
Cushion Sizes	6406-09-A	16" x 16" x 2"	6 x 6	6406-10-A	16" x 16" x 3"	8 x 8
	6406-09-B	18" x 16" x 2"	7 x 6	6406-10-B	18" x 16" x 3"	9 x 8
	6406-09-C	18" x 18" x 2"	7 x 7	6406-10-C	18" x 18" x 3"	9 x 9
				6406-10-D	20" x 16" x 3"	10 x 8
				6406-10-E	20" x 18" x 3"	10 x 9
				6406-10-F	22" x 20" x 3"	11 x 10
Pressure Mapping*	 <p>Average Pressure: 34.6 mmHG</p>			 <p>Average Pressure: 23.0 mmHG</p>		
	Cell Height	2"		3"		
Cushion Material	PVC			Neoprene		
Maximum Weight Capacity	250 lbs.			300 lbs. (Models A-E) 550 lbs. (Model F)		
Cover Available	Yes			Yes		
Cover Washable	Yes			Yes		
Air Pump Standard	Yes			Yes		
Multiple Patient Use	Yes			Yes		
HCPCS Codes	E2622			E2622 & E2623		
Warranty	Eighteen month warranty on manufacturer's defect from the original date of purchase and six month limited warranty for cushion cover					

* Pressure mapping performed with X-Sensor Pressure Mapping equipment. Capillaries close at 32-34 mmHg.