

## Nisus Portfolio Accessories

Cork Medical's wide array of accessories is designed to work in conjunction with the entire portfolio of Nisus Negative Pressure Wound Therapy devices. Each accessory has a simple application process allowing healthcare professionals to choose from multiple application methods designed to promote positive outcomes without sacrificing patient comfort. Take control of the patient's full healing experience.



### Y-Connector

Part Number: **CWY-100**

Cork Medical's Y-Connector is designed to simultaneously treat multiple wounds on the same patient. The multiple wounds should be the same etiology requiring the same pressure settings. The Y-Connector is indicated for use with:

- ▶ All other Nisus NPWT wound dressings



### 24" Extension Tubing

Part Number: **CSP-100-ES-EXT**

Cork Medical's 24" Extension Tubing is used to add length to port pad tubing and will aid in increasing mobility for patients with a wound located on distal portions of the body. The 24" Extension Tubing is indicated for use on:

- ▶ Single wounds where the patient needs longer tubing to increase mobility



### Bridge Foam

Part Number: **NIS-BRDG-A-16**

Cork Medical's Bridge Foam bridges posterior wounds to the lateral or anterior surface, minimizing drainage/port pad tubing-related pressure to surrounding tissue. The 21.75" (L) x 1.96" (W) Cork Medical Bridge Foam is indicated for use when:

- ▶ Bridging the port pad away from any wound where pressure might be applied to the port pad



### Cleanse Therapy Kit

Part Number: **CLEANSE-KIT-A-18**

Cork Medical's NPWT Cleanse Therapy Kit is designed to irrigate the wound while simultaneously providing negative pressure wound therapy. The Cleanse Therapy Kit comes with a port pad and connecting syringe to aid and promote the removal of excess exudates, infectious material, and tissue debris. The kit is ideally used for:

- ▶ Injecting liquids into the wound site for dwelling



## Nisus Portfolio Accessories



### Nisus Canisters

The Nisus Canisters are used in conjunction with the Nisus device and wound dressing kit. The uniquely designed canister collects and stores wound drainage, and its airtight seal allows it to sustain high pressure under any variable condition. This allows patients to maintain flexibility and portability. Nisus Canisters can hold up to two different volume sizes:

- ▶ 250mL for small wounds Part Number: **CPC-250**
- ▶ 500mL for large wounds Part Number: **CPC-500**



### Wound Dressing Kits

The Nisus Wound Dressing Kit's user-friendly design allows for a simplistic application method. Kits come in three different sizes, small, medium, and large, to easily dress wounds for negative pressure wound therapy. Equipped with a hermetic seal, each kit allows for patient mobility as to not sacrifice comfort. Each kit includes:

**Small Kit** Part Number: **WDK-250**

- ▶ A standard size piece of foam
- ▶ One piece of drape
- ▶ A port pad

**Medium kit** Part Number: **WDK-300**

- ▶ A standard size piece of foam
- ▶ Two pieces of drape
- ▶ A port pad

**Large Kit** Part Number: **WDK-400**

- ▶ A large size piece of foam
- ▶ Two pieces of drape
- ▶ A port pad

