

NISUS[®] touch

NEGATIVE PRESSURE WOUND THERAPY PUMP

CORK
MEDICAL

A Unique Combination of Power and Advanced Features

Technology meets performance with the Nisus Touch's intuitive software built for easy usability and maximum healing.



The Cork Medical Nisus Touch Negative Pressure Wound Therapy Pump is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.

NISUS TOUCH NPWT PUMP FEATURES

The Nisus Touch Negative Pressure Wound Therapy Pump is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- ▶ Recommended for patients with chronic, acute, traumatic, sub-acute, dehisced wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, venous ulcers, flaps, and grafts
- ▶ Powerful flow rate of ~8-10 L/min provides the strength needed for healing large and difficult wounds
- ▶ FDA approved for use in hospitals, long term care facilities, and home-care settings
- ▶ Light-weight and portable, with up to 24hrs of battery life
- ▶ User friendly interactive interface with easy to read and touch enabled LCD screen
- ▶ Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage
- ▶ Multiple unique features for clinical enhancement of negative pressure wound therapy such as: patient runtime, alarm log, lock restriction, Variable Intermittent and Continuous Therapy Modes



FDA Approved For Use In Hospitals, Long-Term Care Facilities, And Home-Care Settings



User Friendly Interactive Interface With Easy To Read And Touch Enabled LCD Screen



Step-By-Step Troubleshooting Guide For Common Alarms



Convenient IV Pole Accessibility



Unyielding Flow Rate Power



Designed For Better Clinical Results In Acute And Long-Term Care Facilities



NISUS[®] touch

NEGATIVE PRESSURE WOUND THERAPY PUMP



COMPONENT SPECS

Comprised of three essential components that actively work together to promote wound healing, the Nisus Touch NPWT System includes the Nisus Touch NPWT Pump, Nisus NPWT Canisters, and Cork NPWT Wound Dressing Kits. Please use under the advisement of a physician and carefully read all therapy and device instructions and safety information prior to use.

Use only Cork Medical components with the Nisus Touch Negative Pressure Wound Therapy Pump.

- ▶ Battery charger and carrying bag included

UNIVERSAL NISUS PLATFORM & ACCESSORIES



Nisus Touch NPWT Pump
CMPP-100



Nisus NPWT Canisters
CPC-250 (250ML)
CPC-500 (500ML)



NPWT Wound Dressing Kit
WDK-250 (SMALL)
WDK-300 (MEDIUM)
WDK-400 (LARGE)



NPWT Cleanse Therapy Kit
CLEANSE-KIT-A-18



Bridge Foam
NIS-BRDG-A-16



Optifoam White Foam
FOAM-WHITE-A-16



24" Extension Tubing
CSP-100-ES-EXT



Y-Connector
CWY-100

