

### NISUS<sup>®</sup> NEGATIVE PRESSURE WOUND THERAPY SYSTEM



#### Pump

1. Part # CMPP-100
2. Flow rate of 8-10 LPM (Powerful)
3. Weight of 2 pounds (Portable)
4. Used from the hospital to the home (Continuum of Care)
5. Pressure range of 40 mmHg to 200 mmHg
6. HCPCS code: E2402
7. Advanced troubleshooting menu
8. Icon based menu
9. Sleep mode
10. Up to 48 hour battery life
11. Intermittent and Variable Intermittent mode of therapy

#### Dressing Kits

1. Four sizes available:  
Part # WDK-250: Small Wound Kit, 0.75"x5"x8" Foam, 10"x12" Transparent Drape (1), Port Pad, Ruler  
Part # WDK-300: Medium Wound Kit, 0.75"x5"x8" Foam, 10"x12" Transparent Drape (2), Port Pad, Ruler  
Part # WDK-300T: Thick Wound Kit, 1.25"x5"x8" Foam, 10"x12" Transparent Drape (2), Port Pad, Ruler  
Part # WDK-400: Large Wound Kit, 0.75"x6"x10" Foam, 10"x12" Transparent Drape (2), Port Pad, Ruler
2. HCPCS code: A6550
3. 30 pores per inch wound foam
4. 400-600 micron opening on foam pores
5. Acyclic adhesion on the drape, that is solvent free
6. Soft silicon tubing
7. Port pad designed of soft silicone
8. 1,2,3 application of the drape (Industry standard)

#### Canisters

1. Two sizes available:  
Part # CPC-250: 250-mL Canister Part # CPC-500: 500-mL Canister
2. HCPCS code: A7000
3. Unique one-handed canister attachment
4. 30" drainage with luer lock connection to the port pad tubing