



NEGATIVE PRESSURE WOUND THERAPY SYSTEM

CORK
MEDICAL

Effortless Negative Pressure Wound Therapy Application

Simplify the patient experience and reduce costs without sacrificing performance with the Nisus One.



Cork Medical Nisus One Negative Pressure Wound Therapy System is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris. The single button design makes applying negative pressure wound therapy easier for patients and physicians without compromising the power and performance demanded for optimal wound healing results

NISUS ONE NPWT PUMP FEATURES

The Nisus One Pump is indicated for use in patients who would benefit from negative pressure wound therapy, particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- ▶ Recommended for patients with chronic, acute, traumatic, sub-acute and dehisced wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, and venous ulcers, flaps, and grafts
- ▶ FDA approved for use in hospitals, long term care facilities, and home-care settings
- ▶ Single button operation provides both patients and physicians with an easy to use, intuitive design
- ▶ Variable Intermittent mode allows users to apply high and low pressure settings at 125 mmhg over 40 mmhg
- ▶ Continuous Therapy mode allows users to set one level of pressure at 125 mmhg
- ▶ Light-weight and portable, weighing less than 2 lbs. with up to 24 hours battery life



Single Button Operation Provides Both Patients And Physicians With An Easy To Use, Intuitive Design



Lightweight And Portable Allowing For Uninterrupted Patient Mobility



FDA Approved For Use In Hospitals, Long-Term Care Facilities, And Home-Care Settings



Multiple Therapy Modes



Leading Affordable Solution



Simplifies The User Experience





NEGATIVE PRESSURE WOUND THERAPY SYSTEM



COMPONENT SPECS

The Nisus One Pump is part of a complete wound care system. Two essential components that actively work together to promote wound healing are the Nisus Canisters and the Cork NPWT Wound Dressing Kit.

Please use under the advisement of a physician. Carefully read all therapy and device instructions and safety information prior to use. Use only Cork Medical components with the Nisus One Negative Pressure Wound Therapy Pump.

- ▶ Battery charger and carrying bag included

UNIVERSAL NISUS PLATFORM & ACCESSORIES



Nisus One NPWT Pump
CMPP-100



Nisus NPWT Canisters
CPC-250 (250ML)
CPC-500 (500ML)



NPWT Wound Dressing Kit
WDK-250 (SMALL)
WDK-300 (MEDIUM)
WDK-400 (LARGE)



NPWT Cleanse Therapy Kit
CLEANSE-KIT-A-18



Bridge Foam
NIS-BRDG-A-16



Optifoam White Foam
FOAM-WHITE-A-16



24" Extension Tubing
CSP-100-ES-EXT



Y-Connector
CWY-100

