

Your First Choice for Wound Healing.

Cork Medical is changing the perception of negative pressure wound therapy (NPWT) with its Versa NPWT device. An ideal solution for healing wounds, the Versa is creating superior patient outcomes by offering NPWT as an affordable and user-friendly solution.



Cork Medical's Versa NPWT device is designed to provide top clinical solutions by maximizing the patient's healing outcome while cutting the cost of managing NPWT.

THE VERSA SYSTEM FEATURES

Cork Medical's Versa Negative Pressure Wound Therapy System is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- ▶ Continuous and intermittent modes to customize patient's therapy
- ▶ Accurate pressure tolerance at wound base
- ▶ Easy to release 100ml canister
- ▶ Standard wound dressing kit for clinical convenience
- ▶ On screen alarm troubleshooting
- ▶ One week of battery life with three easily replaceable AA lithium batteries
- ▶ Lightweight and portable, weighing less than one pound
- ▶ FDA approved for use in long term care facilities, acute care facilities, and home care



The Versa's low entry cost allows facilities to use the device as a one-time use or up to multiple uses, reducing the need for inventory management and service or repairs.



A compact and simple design allows clinical staff to easily administer negative pressure.



Its low cost and superior healing abilities allow facilities to utilize the healing power of NPWT as their first-choice treatment option.



VERSA



COMPONENT SPECS

Comprised of three essential components that actively work together to promote wound healing, the Versa NPWT System includes the Versa NPWT Pump, the Versa NPWT Canisters, and Cork Medical's NPWT Wound Dressing Kit. Please use under the advisement of a physician. Carefully read all therapy and device instructions and safety information prior to use.

Use only Cork Medical components with the Versa Negative Pressure Wound Therapy System.

- ▶ Batteries and carrying case included

VERSA NPWT SYSTEM ACCESSORIES



Versa 100ml Canisters
CPC-100 (100ML)



AA Lithium Batteries
VRSA-POWER-A-21



NPWT Wound Dressing Kit
WDK-250 (SMALL)
WDK-300 (MEDIUM)
WDK-400 (LARGE)



Versa Carrying Case
CWN-100

