

## Forever Plus Cushion

Nineteen individual air cells provide a unique design which keeps the lower layer of air cells constantly inflated while alternating cells inflate and deflate the upper layer.



### KEY BENEFITS OF THE FOREVER PLUS CUSHION

The Forever Plus Cushion Series offers a therapeutic cushioning design, constructed of independent flexible interconnected air cells, dispersing pressure evenly to increase blood flow and assist in the treatment and prevention of pressure ulcers. Flexible 3-inch interconnected independent air cells allow for variable patient immersion into the cushion, providing consistent low pressure and an increased blood supply necessary for healthy tissue maintenance and treatment.

- ▶ The water resistant, breathable cover provides a cool, comfortable seating environment
- ▶ Multi-cellular, air-filled design allows for adjustments and a personalized fit
- ▶ Facilitates blood flow by fitting and conforming to an individual's shape

### AVAILABLE IN 6 SIZES

Model	Height	Dimensions	Cells	Material	Max. Wt
6406-10-A	3"	16" x 16"	8 x 8	Neoprene	300 lbs.
6406-10-B	3"	18" x 16"	9 x 8	Neoprene	300 lbs.
6406-10-C	3"	18" x 18"	9 x 9	Neoprene	300 lbs.
6406-10-D	3"	20" x 16"	10 x 8	Neoprene	300 lbs.
6406-10-E	3"	20" x 18"	10 x 9	Neoprene	300 lbs.
6406-10-F	3"	22" x 20"	11 x 10	Neoprene	550 lbs.

HCPCS Code: **E2622**

### COMPONENT SPECS

- Machine washable cover has non-slip coating on bottom to prevent sliding
- Includes air injection hose, hand-held pump, inflation valve, and cushion cover
- Weight capacity: 300 lbs. (Models 6406-10-A through E), 550 lbs. (Model 6406-10-F)
- Medicare HCPCS Codes: E2622 (Models 6406-10-A through E), E2623 (Model 6406-10-F)

