

NISUS[®] touch NEGATIVE PRESSURE WOUND THERAPY PUMP



The Cork Medical Products Nisus Touch Negative Pressure Wound Therapy Pump is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.

NISUS TOUCH NPWT PUMP FEATURES

The Nisus Touch Negative Pressure Wound Therapy Pump is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- Recommended for patients with chronic, acute, traumatic, sub-acute, dehisced wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, venous ulcers, flaps, and grafts
- Powerful flow rate of ~8-10 L/min provides the strength needed for healing large and difficult wounds
- FDA approved for use in hospitals, long term care facilities, and home-care settings
- Light-weight and portable, with up to 24hrs of battery life
- User friendly interactive interface with easy to read and touch enabled LCD screen
- Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage
- Multiple unique features for clinical enhancement of negative pressure wound therapy such as: patient runtime, alarm log, lock restriction, Variable Intermittent and Continuous Therapy Modes

COMPONENT SPECS

Comprised of three essential components that actively work together to promote wound healing, the Nisus Touch NPWT System includes the Nisus Touch NPWT Pump, Nisus NPWT Canisters, and Cork NPWT Wound Dressing Kits.

Please use under the advisement of a physician and carefully read all therapy and device instructions and safety information prior to use.

Use only Cork Medical Products components with the Nisus Touch Negative Pressure Wound Therapy Pump.

Model No.	Component
TCMPP-100	Nisus Touch NPWT Pump
CPC-250	Nisus NPWT 250mL Canister
CPC-500	Nisus NPWT 500mL Canister
WDK-250	NPWT Wound Dressing Kit (Small)
WDK-300	NPWT Wound Dressing Kit (Medium)
WDK-400	NPWT Wound Dressing Kit (Large)
CLEANSE-KIT-A-18	NPWT Cleanse Therapy Kit
NIS-BRDG-A-16	Bridge Foam
FOAM-WHITE-A-16	HydroLite White Foam
CSP-100-ES-EXT	24" Extension Tubing
CWY-100	Y-Connector

- Battery charger and carrying bag included
- Medicare HCPCS Code: E2402

**CORK
MEDICAL**


1-866-551-2580

Revision 2.0 2021_1006

NISUS® touch

NEGATIVE PRESSURE WOUND THERAPY PUMP

User friendly interface with easy to read and touch enabled LCD screen



Latch design provides streamlined, one-handed canister release and removal



The Cork Medical Products Nisus Touch Negative Pressure Wound Therapy Pump is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.

DEVICE FEATURES

The Cork Nisus Touch Negative Pressure Wound Therapy Pump is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- Recommended for patients with chronic, acute, traumatic, sub-acute, dehisced wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, venous ulcers, flaps, and grafts
- Powerful flow rate of ~8-10 L/min provides the strength needed for healing large and difficult wounds
- FDA approved for use in hospitals, long term care facilities, and home-care settings
- Light-weight and portable, weighing less than 2 lbs. with up to 24 hours battery life
- Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage
- Multiple features for clinical enhancement such as: patient runtime, alarm log, lock restriction, Variable Intermittent and Continuous Therapy Modes

DEVICE SPECS

The Nisus Touch Negative Pressure Wound Therapy Pump is part of a complete wound care v. Two essential components that actively work together to promote wound healing are the Nisus NPWT Canisters and the Cork NPWT Wound Dressing Kit.

Please use under the advisement of a physician and carefully read all therapy, device instructions and safety information prior to use.

Use only Cork Medical Products components with the Nisus Touch Negative Pressure Wound Therapy Pump.

- Battery charger and carrying bag included
- Medicare HCPCS Code: E2402

Model No.

TCMPP-100
CPC-250
CPC-500
WDK-250
WDK-300
WDK-400
CLEANSE-KIT-A-18
NIS-BRDG-A-16
FOAM-WHITE-A-16
CSP-100-ES-EXT
CWY-100

Component

Nisus NPWT Pump
Nisus NPWT 250mL Canister
Nisus NPWT 500mL Canister
NPWT Wound Dressing Kit (Small)
NPWT Wound Dressing Kit (Medium)
NPWT Wound Dressing Kit (Large)
NPWT Cleanse Therapy Kit
Bridge Foam
HydroLite White Foam
24" Extension Tubing
Y-Connector



1-866-551-2580

Revision 2.0 2021_1006