

The True Balance Air is a highly effective, non-powered pressure redistribution mattress with combined high density foam and self-adjusting air cells, providing optimum support for the treatment and prevention of pressure ulcers.



Posterior high-density memory foam layer prevents bottoming-out from occurring

Eight self-adjusting air cells for maximum pressure redistribution

### FEATURES OF THE TRUE BALANCE AIR SERIES

Eight horizontal, interconnected air chambers automatically adjust to a patient's position for optimum pressure redistribution.

- Anti-Microbial and hypo-allergenic replaceable cover is easy to clean
- Zipper-flap closure for easy access to mattress
- Posterior high density foam layer prevents bottoming-out from occurring
- Interconnected air cells for even, quiet, non-powered pressure redistribution
- Sloped heel design for optimal heel pressure reduction
- Standard fire barrier for regulatory compliance

### SPECS — STANDARD, PART NUMBER H12088

- Size of mattress: 80" (L) x 36" (W) x 7" (H)
- 8 air cells: 25.5" (L) x 6" (W) x 3.5" (H)
- Maximum weight capacity: 350 lbs.
- Material: Urethane cover and Tri-Laminate vinyl bottom cover; High density, anti-microbial foam
- 1-Year Warranty
- Medicare HCPCS Code: E0373

### SPECS — HD, PART NUMBER 6406-12

- Size of mattress: 80" (L) x 42" (W) x 11" (H)
- 8 air cells: 33" (L) x 6" (W) x 3.5" (H)
- Maximum weight capacity: 600 lbs.
- Material: Urethane cover and Tri-Laminate vinyl bottom cover; High density, anti-microbial foam
- 1-Year Warranty
- Medicare HCPCS Code: E0373