

## Elite Balance Air Series

The Elite Balance Air Series is an alternating pressure mattress replacement system for use in the treatment and prevention of pressure ulcers. It provides patients a comfortable and therapeutic support surface to both enhance healing and prevent skin breakdown.



Standard and HD Mattresses are available with or without bolsters

Replacement covers available

## Features of the Elite Balance Air

Nineteen individual air cells provide a unique design which keeps the lower layer of air cells constantly inflated while alternating cells inflate and deflate the upper layer.

- Digital pump with light-up display allows for effortless user experience and control
- Pressure control settings allow patient to set their personalized comfort setting depending on weight
- Auto-Detection setting automatically detects the patient weight and sets the appropriate pressure output for optimum patient comfort and ease
- Pump features audible and visual low pressure, power failure and machine malfunction alarms
- Split cell-on-cell design allows the patient to remain in an upright position and keeps mattress inflated up to 24 hours if loss of power should occur
- In the event of a medical emergency, mattress features a rapid deflation CPR valve

## Specs

- 19 air cells
- Cell Material: Nylon TPU
- Cover Material: Stretched Polyester with Polyurethane
- 1-Year Warranty
- Medicare HCPCS Code: E0277

## Available in 3 Sizes

Model Number 6406-05 or 6406-05-Bolster 6406-08 or 6406-08-Bolster 6406-11 Model Name EBA Standard EBA HD EBA MAX Dimensions 78" (L) × 35" (W) × 8" (H) 78" (L) × 42" (W) × 8" (H) 78" (L) × 48" (W) × 8" (H)

Max. Weigh
350 lbs.
550 lbs.
750 lbs.

