

The Comfort Select Gel Overlay, Group 1 Support Surface provides superior pressure redistribution and optimal patient comfort.



Unique poly-gel core design lays easily on top of existing mattress

15 gel bladders, each with 3 individual channels, 2" thick allowing 75% surface area coverage

High-density, flame-retardant polyurethane foam casing surrounds the gel bladders

FEATURES OF THE COMFORT SELECT GEL

Fifteen gel-filled bladders, each with three individual channels, cover 75% of the surface area for maximum client comfort, while providing optimum pressure redistribution.

- Convenient one-piece construction requires no assembly
- Reduce shearing injuries through easy client transfer with an anti-shear top cover
- Anti-bacterial, waterproof top cover is flexible and flame retardant
- 4 elastic corner straps ensure the overlay is secured to existing mattress

SPECS

- Size of overlay: 76" (L) x 34" (W) x 3.5" (H)
- 15 gel bladder cells (each with 3 individual channels): 65" (L) x 30" (W) x 2" (H)
- Maximum weight capacity: 275 lbs.
- Cover Material: Soft nylon urethane-coated, non-vapor material; vinyl bottom cover
- Overlay Material: High-density, polyurethane foam casing
- 1-Year Warranty
- Medicare HCPCS Code: E0185